

Voice Commands: 09-14-2021

OrCam has 4 Voice Command Groups:

- 1- Smart Reading
- 2- Orientation Commands
- 3- Information (Tell) Commands
- 4- Settings Commands

*You can also ask OrCam by using **“Hey OrCam”** to list voice commands and say which group you want

SMART READING	SETTINGS	INFORMATION (TELL)	ORIENTATION
Smart reading	Speak Louder / Softer	Tell the time	What's in front of me
Read everything	Read Faster / Slower	Tell the date	Repeat
Find	Enable / Disable the Reading Navigation	Tell the battery status	Exit
Start from	Set the face recognition to auto	Tell the device serial ID number	
Jump too	Set the face recognition to manual	Tell the Version	
Read the article about	Disable the face recognition		
Read the dates	Enable / Disable the Auto page recognition		
Read the phone numbers	Enter the user menu		
Read the amounts	Change the voice		
Previous	Change the auto suspend time		
Next	Change the face repletion time		
Repeat	Start the face learning		
Read the headlines	Start barcode learning		
Read article number 1-10	Start product learning		
Read the first article / second article...	Enable /Disable the auto flash		
Exit	Connect to the Bluetooth audio device		
Restart	Check the internet connectivity		
Volume up / down			
Speak louder /softer			
Reading speed up / down			
Read faster / slower			
Help			